

Test Order

Shiga Toxin Producing Escherichia coli (STEC) Toxin Testing

Shiga toxins 1 and 2 Immunoassay
Immunochromatographic Rapid Test
Pure isolates, previously identified as Escherichia coli
Transport on nutrient agar slant, packaged and sealed according to
regulations.
Refrigerate until transported to the laboratory. Transport at room
temperature.
N/A
N/A
N/A
No
7 business days (M-F)
http://sfcdcp.org/document.html?id=1036

Updated 2/15/2017